



Model: E817U

mailbox: admin@yanranyl.com

Company: Guangzhou Yanran Medical Equipment Technology Co., Ltd.

Address: Room 401, Building G, No. 6 Helong Yiheng Road,
Baiyun District, Guangzhou

Company website: www.yanranyl.com

Ion detox foot spa Operation Manual

Please read through this manual before usage.
By using you agree to have read the manual.

Basic parameters

product name: Ion detox foot spa

Input: AC110-220V~50Hz

Output: DC12V

Output power: 24W



control unit

Matching set



power supply

array

spoon

- foot tub--authentication by SGS
- control unit
- array
- power supply
- spoon

Basic parameters

1. Place a PE bag in the foot tub and fill with warm Water .
2. Add 0.1g sodium salt into the water.(spoon supplied in pack)
3. FiX the control unit to the foot tub.(Please see the picture below)
4. Place the array in the center of the foot tub and connect to the control unit .
5. Connect the power supply to the control unit and plug switch into the main machine.
6. Put feet into foot tub,push button on control panel.
7. The series of lights on the control panel illuminate during the work progresses,Each light illuminate for 5 minutes,the treatment will last 30 minutes. (The lights that are already on will stay on until the end of the work)
8. A buzzer will sound intermittently at the end of 30minutes and all the lights will flash four times to Indicate work has finished.



Indicator light
(lighting state)




control unit

switch



Color reference



	Yellow-green	Detoxifying from the kidney, bladder urinary tract, female prostate
	Orange	Detoxifying from joints
	Brown	Detoxifying liver, tobacco, cellular debris
	Black	Detoxifying liver, gallbladder
	Dark Green	Detoxifying gallbladder
	White Foam	Mucus from lymph
	White Speckles	Mostly yeast
	Black Speckles	Heavy metals
	Red Speckles	Blood clot material

Explanation

1. During the 30 minutes the foot spa will begin to rebalance the positive and negative ions in your body. It is important to drink at least 500ml of water during treatment.
2. You will notice the water will change colour through treatment, this is perfectly normal. You may also see deposits floating on top of the water. Again this is normal as one of the benefits of a foot spa treatment is Lymphatic drainage.
3. The array will usually last between 30 and 50 treatments. When the array stops working there will be no color change to the water.

Caution

Please note as following:

1. If low blood pressure, please eat something first. The machine may decrease blood sugar.
2. The pregnant or lactation women mustn't use it, because it maybe dredge out the toxins of the tissue and grease and block producing milk nutrition that is necessary to the infant body. Although it not be proved with bad impaction, but for safety, we suggest not to use it.
3. After use, please make up nutrition absorption such as vitamin, mineral, fatty acid and etc to promote the body metabolism and eject the toxin by itself.
4. To avoid any effects of electromagnetic interference do not use the appliance near a mobile phone or any other powerful radio interference equipment.

Questions and answers

Will I feel any sensation during the treatment?

This depends some people do experience tingling sensations, these can be in your feet, legs or arms.

Will I suffer any side effects from the treatment?

Generally no. However, some people have experienced mild headaches, or mild flu like symptoms, after treatment. Drinking plenty of water during the treatment cycle will reduce these symptoms or totally prevent them.

Do I need to do anything during the treatment?

No. Just sit down and relax for the duration of the Treatment.

How many treatments should I have?

Ideally two per week. However, there must be at least 48 hours between treatments.

Is there anyone who should not have treatments?

Anyone who suffers from epilepsy, has a pacemaker fitted, has had organ implants, if you are pregnant, under the age of 10, has any open wounds or abrasions on their feet should not have treatments.

Why does there appear to be deposits floating on top of the Water?

Again, this can be as a result of the reaction between the array and saline solution. Or, it can be deposits drained from the lymph glands and or pores.

How long will the array last?

Usually between 30 and 50 times. When the array 'expires' it will simply stop working.

Why does the Water change color?

There will be a reaction between the array and the saline solution, which will give a light orange/yellow tint to the water. This is a natural occurrence and is the sodium being drawn out of the salt by the current.

Questions and answers

We live in a hard water area will this have any effect on the treatment?

Yes, you should reduce the amount of salt you add to the water.

If the machine shut off automatically, what should I do?

First, please check the connection and make sure that it connects well; Second, replace the water, as the salt is too much, the machine will stop working automatically.

Who should not use the Ion-Cleanse?

The Ion-Cleanse is safe to use by most people, except maybe in the case of

- During pregnancy (or suspected pregnancy)
- People with a pacemaker
- People with implanted organs or steel plates or joints
- People with epilepsy
- People undergoing any form of radiotherapy or chemotherapy
- People with type I Diabetes
- People with open wounds on their feet
- If you are unsure please check with your G.P or Specialist.